

Contec® Sterile 70% Isopropanol

USP <797> Compliant Sterile Alcohol

Contec Sterile 70% Isopropanol contains 70% by volume USP Grade Isopropanol (IPA) and 30% USP Grade Purified Water. Filtered to 0.2µ and induction sealed to assure sterility, the bottles are double-bagged and gamma-irradiated. Each bottle has a lot code and expiration date for easy record keeping. Product expires 5 years after date of manufacture. Validated per the Association for the Advancement of Medical Instrumentation (AAMI) 11137:2006 Guidelines to a 10⁻⁶ Sterility Assurance Level (SAL). *Note: Sterility validation is performed after trigger is installed.*

A convenient trigger sprayer with settings for stream or mist spraying is included with 16 and 32 oz bottles. The stream setting is ideal for targeting specific hard to reach areas, while the mist setting is used to cover larger areas. Each leak-proof bottle is induction sealed, capped and shipped with an uninstalled trigger sprayer to eliminate the occurrence of leakage. Instructions for sterile installation of the trigger sprayer are available.



Applications:

- Wipe down of isolators, cabinets, laminar flow hoods and surfaces
- Use on gloved hands in sterile suites
- Wipe down of products and equipment for pass through to controlled environments
- Ideally suited for use with sterile cleanroom wipes

Appearance	Clear, colorless liquid
Odor	Characteristic alcohol
Water Quality	Meets the requirements of USP Grade Deionized Water
% IPA Contents (+/-3%)	67.9% - 72.1%
Specific Gravity (20°C)	0.872 - 0.883

FEATURE	BENEFIT
Contains 70% USP-grade Isopropanol	<ul style="list-style-type: none"> • Compliant with USP <797> • Evaporates and leaves extremely low residues • Eliminates mixing, filtration, processing and QC expense of in-house production of sterile Isopropanol water blends
Filtered through 0.2 micron in-line filter	<ul style="list-style-type: none"> • Eliminates the need for in-house filtration • Meets USP Grade Purified Water
Gamma-irradiated to a sterility assurance level (SAL) of 10 ⁻⁶	<ul style="list-style-type: none"> • Compliant with USP <797> • Prefiltered, packaged and sterilized for easy introduction into critical environments • Certificate of processing and sterility with each case
Induction sealed bottles	<ul style="list-style-type: none"> • Reduces opportunity for leakage during sterilization and shipment
Adjustable trigger spray	<ul style="list-style-type: none"> • Select spray or jet option
Plastic bottle	<ul style="list-style-type: none"> • Recyclable, environmentally friendly, no costly disposal
Double bagged	<ul style="list-style-type: none"> • To preserve product integrity
Lot numbered bottles	<ul style="list-style-type: none"> • Assures complete product traceability

NOT FOR RETAIL SALE. SAFETY DATA SHEET AVAILABLE FOR PROFESSIONAL USER UPON ON WWW.CONTECINC.COM.

Part No.	Description	Quantity
SB1670301R	Sterile 70% Isopropanol, 16 oz. (0.47L) w/trigger sprayer	12 bottles/case
SB3270301R	Sterile 70% Isopropanol, 32 oz. (0.95L) w/trigger sprayer	12 bottles/case
SB12870301R	Sterile 70% Isopropanol, 1 gallon (3.8L) bottle	4 bottles/case



P.O. Box 530 • Spartanburg, SC 29304 USA
 Tel: +1-864-503-8333
 Toll Free (US, Canada, PR): 1-800-289-5762
 Europe: +33 297437690
 www.contecinc.com • disinfect@contecinc.com

* Also available nonsterile. See website for product offering.